

# Futera Dots

Non-Invasive Fractional RF

## 膠原復活槍 重拾0毛孔

非侵入性 | 無副作用 | 快速復原

透過非侵入性多極分段式射頻(UFMR)·  
釋放熱力「氣化」肌膚表皮而獲得重新·  
在修復過程會促進更新及刺激膠原蛋白·  
讓皮膚緊緻·重拾零毛孔示人。



UFMR (Ultra Fractional Multi-Polarity Radiofrequency)


### 多極分段式射頻技術

以1 MHz 聚焦分段式射頻·釋放熱力達到  
55~60°C·瞬間「氣化」肌膚表皮·能量  
有效輸送真皮層·加熱深度小於400微米·  
確保以最小熱點達輕微創傷·令肌膚自行  
更生·促進及刺激膠原蛋白增生·在復原  
過程中會拉扯表皮·收緊鬆弛皮膚和收細  
毛孔·達到【先死而後生】!

傳統  
分段激光      UFMR  
多極分段式射頻



多極64Pin細小機頭·接觸式傳送·全面照顧鼻翼、下巴、凹陷、皺紋



電極接觸皮膚·因溫度產生熱損傷·點與點之間分  
隔保留原有皮膚·有效迅速康復而創傷性少

技術	多極分段式射頻	微針點陣射頻	分段式二氧化碳激光
運作	正 / 負電極組成，接觸表皮 	通過微針電極刺入皮膚 	整層表皮氧化 (磨去) 
創傷	最小熱點，輕微創傷	會有針孔，中度創傷	磨蝕性，極度創傷
復原期	1-2天	1星期	1-2個月
即時反應	微紅，溫熱，淺層點狀結痂	微紅，溫熱，點狀結痂	維持1-2天紅腫，灼熱感，無數微小的點狀傷口
安全	不用麻醉，不會出血	可麻醉，輕微滲血	須麻醉
皮膚種類	任何皮膚	任何皮膚	不能使用在深色皮膚

### Whitening



Before

After


### Tightening Pores



Before

After


### Wrinkle Improvement



Before

After

### Skin tone Improvement



Before

After

## Specification

Country	Korea
Product Type	Fractional RF
Frequency Output	1MHz
Tip type	64Pin
Output Range	80W (Rated Load 200 Ω)
Repeat time	0.8~1.4 Sec.
Exposure time	5~200ms
Treatment area	11x11 mm
Electrical Req.	AC110~220V,50/60Hz
Dimension(WxDxH)	340x350x280(mm)
Weight	4.5kg
User Interface	8" Touch LCD
Power Consumption	180VA



### 4次療程，3-4週一次

- 撫平凹凸洞，平滑肌膚
- 減少黑頭粉刺，收細毛孔
- 緊緻鬆弛表皮
- 提升輪廓，增加彈性

- 減淡皺紋及幼紋
- 改善痤瘡，疤痕，暗瘡印
- 均勻膚色，增加光澤
- 改善妊娠紋
- 療程後隔日可上妝